

*Internal use-only*  
**Work Order #** \_\_\_\_\_  
 D.R. \_\_\_\_/\_\_\_\_/\_\_\_\_  
 D.C. \_\_\_\_/\_\_\_\_/\_\_\_\_  
 D.S. \_\_\_\_/\_\_\_\_/\_\_\_\_  
 Scan #: \_\_\_\_\_

## LAB CLIENT INFO:

### \*Required Information In Shaded Area:

Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 \_\_\_\_\_  
 Phone #: \_\_\_\_\_

Patient Name: \_\_\_\_\_  
 Shoe size: \_\_\_\_\_ Weight: \_\_\_\_\_  
 Footwear type: \_\_\_\_\_ Male / Female (circle)  
 Activities: \_\_\_\_\_

*\*Please Check Off Only ONE Orthotic Per Order Form*

### Standard Orthotics

- Day-to-Day Orthotic (Casual shoes)**  
Basic day-to-day orthotic for most standard footwear.
- Children's Orthotic**  
Basic day-to-day orthotic for most standard children's footwear.

### Dress Orthotics

- Women's / Men's Dress Orthotic (Flat dress shoes)**  
Slim design with shallow heel cup for low flat lace up or slip on dress shoes.
- Women's Heels (Cobra Design)**  
Designed for heel shoes of heel height from 1/2" – 1 1/2".

### Sport Orthotics

- Regular Sport Orthotic (Athletic shoes)**  
Designed to work with most standard athletic shoes.
- Low Profile Sport Orthotic (Ski/Skate/Cleat/Cycling)**  
Designed to work with footwear that requires less volume.
- Max Impact Sport Orthotic**  
Designed to work in running shoes and boots for extra cushioning.








### Soft Orthotics

- Soft Shell XPE**  
2 layer design for increased cushioning and flexibility for sensitive feet.
- Diabetic / Arthritic Orthotic**  
Provides maximum cushioning and pressure relief. Requires a shoe with more volume.

### Shell Material and Support Rating

Shell Material:  Polypropylene  Sub-ortholene  Graphite  
 Support Rating:  Flexible  Semi-Flex  Semi-Rigid  Rigid  
 Notes: \_\_\_\_\_

### Additions/Modifications

 Heel Lift: \_\_\_\_\_ mm / inch Attached / Separate L / R  
 Heel Well: \_\_\_\_\_ L / R  
 Metatarsal Pad (circle one): Low Medium High L / R  
Small Medium Large L / R  
 Metatarsal Bar: \_\_\_\_\_ L / R  
 Neuroma Pad: Placement 1 2 3 4 5 L / R  
 Morton's Extension: Shell / Extrinsic L / R  
 Reverse Morton's: Shell / FF Ext / Both L / R




### Top-covers:

- Vinyl  1/16" MC Puff
- Leather  1/8" MC Puff
- Neoprene  X-Static
- Plastazote  UltraHyde

### Cover Notes:

\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

### Top-cover length:

- 3/4 Length 
- Sulcus 
- Full 

### Posting Instructions

Location	Posting	Left	Right
Hindfoot	Internal / External	Var/Val _____°	Var/Val _____°
Forefoot	Internal / External	Var/Val _____°	Var/Val _____°

### Special Instructions

\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

### Lab Use Only

\_\_\_\_\_  
 \_\_\_\_\_